


SAMPLE INFORMATION

Sample Name:	M40701-Suzetrigine-02-1	Acquired By:	System
Sample Type:	Unknown	Sample Set Name:	m40701
Vial:	8	Acq. Method Set:	CD_05_95_15min_40D
Injection #:	1	Processing Method:	LC PQ
Injection Volume:	5.00 ul	Channel Name:	220.0nm
Run Time:	15.0 Minutes	Proc. Chnl. Descr.:	W2996 PDA 220.0 nm (PDA)



	RT	Area	% Area	Height
1	7.249	13363	0.09	3271
2	8.959	14129220	99.84	2454170
3	10.573	9151	0.06	1883