




Agilent Technologies

Sample name: M40714-N-Lactoyl-Phenylalanine-02 Injection date:
Instrument: A012 Location:
Inj. volume: 2.5 Type:
Acq. method: HPLC Sample amount:
Processing method:
Manually modified:



Signal: DAD1B,Sig=210.0,4.0 Ref=off

RT [min]	Area	Height	Area%
5.770	4.14	1.57	0.05
5.936	7555.19	1965.00	99.71
6.392	5.57	2.19	0.07
6.468	5.09	1.97	0.07
7.056	6.79	2.37	0.09