




Agilent Technologies

Sample name: M4138-Deucravacitinib (BMS-986165)-02 Injection date:
Instrument: A010 Location:
Inj. volume: 2 Type:
Acq. method: HPLC Sample amount:
Processing method:
Manually modified:



Signal: DAD1D, Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
4.944	BB	0.17	11.95	3.66	0.12
5.176	MM m	0.10	1.40	0.54	0.01
5.317	BM m	0.16	23.67	7.37	0.24
5.436	MM m	0.08	1.31	0.58	0.01
5.616	BB	0.18	24.98	6.76	0.25
5.873	BM m	0.35	9827.85	2274.70	99.17
6.196	MM m	0.10	2.79	0.52	0.03
6.315	MM m	0.14	1.36	0.29	0.01
6.701	MM m	0.20	5.97	0.82	0.06
6.887	MM m	0.10	1.08	0.36	0.01
7.006	MM m	0.14	1.70	0.47	0.02
7.440	BB	0.13	4.71	1.34	0.05
7.681	MM m	0.12	1.51	0.45	0.02