



Sample name: M4139-Salvianolic-acid-A-02

Injection date:

Instrument: A010

Location:

Inj. volume: 2


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1B, Sig=210.0,4.0 Ref=off

RT [min]	Area	Height	Area%
6.127	12.52	2.55	0.14
6.242	3.70	0.83	0.04
6.484	8983.57	1825.24	98.93
6.810	18.99	5.01	0.21
6.936	62.14	14.36	0.68