




$^1\text{H}$  NMR (400 MHz,  $\text{dms-}d_6$ )  $\delta$  10.06 (s, 1H), 9.78 (s, 1H), 8.66 (s, 1H), 8.47 (s, 1H), 7.79 (s, 1H), 7.57 – 7.44 (m, 2H), 7.28 (td,  $J$  = 8.1, 1.3 Hz, 1H), 6.80 (dt,  $J$  = 15.4, 5.9 Hz, 1H), 6.44 (d,  $J$  = 15.5 Hz, 1H), 3.61 (t,  $J$  = 4.6 Hz, 4H), 3.19 – 3.06 (m, 3H), 2.92 (d,  $J$  = 9.1 Hz, 1H), 2.48 – 2.32 (m, 6H), 2.25 (s, 3H), 1.93 (dt,  $J$  = 8.2, 4.4 Hz, 1H), 1.37 (t,  $J$  = 4.3 Hz, 1H), 1.05 – 1.00 (m, 1H).

