



Agilent Technologies

Sample name: M4240-Shikonin-03

Instrument: A012

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B,Sig=210.0,4.0 Ref=off

RT [min]	Area	Height	Area%
8.616	5.65	1.48	0.04
8.831	8.50	2.27	0.05
9.414	2.36	0.57	0.02
10.531	15521.97	2008.15	99.58
16.721	24.64	4.51	0.16
17.047	25.07	5.88	0.16