




Sample name: M42874-JY-2-01
Instrument: A010
Inj. volume: 1.5
Acq. method: HPLC
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1A, Sig=220,4 Ref=off

RT [min]	Width [min]	Area	Height	Area%
5.731	0.19	8.3	2.6	0.047
6.975	0.24	107.9	27.4	0.606
7.619	0.13	1.7	0.5	0.010
8.302	0.20	4.3	0.9	0.024
8.885	0.22	15.8	1.9	0.089
9.073	0.15	9.6	2.3	0.054
9.161	0.14	27.8	6.4	0.156
9.522	0.41	61.0	8.4	0.343
10.556	0.76	17547.8	2256.5	98.607
11.999	0.16	11.4	2.6	0.064