



Agilent Technologies

Sample name: M4394-Wogonoside-02

Instrument: A010

Inj. volume: 1.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B, Sig=210.0, 4.0 Ref=off

RT [min]	Area	Height	Area%
4.936	4.63	1.73	0.08
6.692	38.80	10.00	0.65
6.817	5957.68	1561.90	99.06
7.402	12.93	4.14	0.21