



DMSO, F1: 400.182, F2: 1.000, SW1: 8621, OF1: 1900.0, PTS1d: 65536, EX: zg30, PW: 12.0 usec, PD: 1.0 sec, NA: 8, LB: 0.0, Sample Name: M4460-Oroxylin-A-01

F1: 400.182	F2: 1.000	SW1: 8621	OF1: 1900.0	PTS1d: 65536
EX: zg30	PW: 12.0 usec	PD: 1.0 sec	NA: 8	LB: 0.0

Sample Name: M4460-Oroxylin-A-01