


SAMPLE INFORMATION

Sample Name:	M4498-Hydroxysafflor-yellow-A-02	Acquired By:	System
Sample Type:	Unknown	Sample Set Name:	M4498
Vial:	87	Acq. Method Set:	CD_240809H_150_AQ
Injection #:	1	Processing Method:	LC PQ
Injection Volume:	10.00 ul	Channel Name:	380.0nm
Run Time:	15.0 Minutes	Proc. Chnl. Descr.:	W2996 PDA 380.0 nm (PDA)



	RT	Area	% Area	Height
1	6.308	2882	0.07	832
2	6.773	3781989	98.35	976613
3	6.947	15494	0.40	6530
4	7.046	38561	1.00	5860
5	7.449	5376	0.14	1574
6	7.520	1333	0.03	517