



Sample name: M45136-PF-07220060-02

Injection date:

Instrument: A012

Location:

Inj. volume: 2


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1E,Sig=240,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
4.440	MM m	0.17	1.09	0.18	0.01
5.056	MM m	0.12	1.24	0.42	0.01
5.506	MM m	0.11	2.32	0.72	0.03
6.008	MM m	0.10	1.55	0.57	0.02
6.154	MM m	0.17	2.74	0.75	0.03
6.378	BB	0.11	25.93	8.90	0.29
6.542	BB	0.37	8857.65	2716.72	99.42
6.839	BB	0.09	3.23	1.21	0.04
7.089	MM m	0.12	1.95	0.57	0.02
7.436	MM m	0.12	2.00	0.60	0.02
8.165	MB m	0.22	9.35	1.40	0.10