



Agilent Technologies

Sample name: M4518-Procyanidin-B2-01

Injection date:

Instrument: 2

Location:

Inj. volume: HPLC


Type:

Acq. method:

Sample amount:

Processing method:

Manually modified:



Signal: DAD1D, Sig=235,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%	Name
4.230	MB m	0.06	5.56	1.29	0.04	
4.471	MB m	0.07	12692.29	2764.43	99.26	
4.695	BB	0.17	89.60	23.00	0.70	
Sum			12787.46			