




Sample name: M4574-Wedelolactone-05
Instrument: A010
Inj. volume: 4
Acq. method: HPLC
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1A,Sig=220.0,4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
6.863		11.1	2.9	0.075
7.150		0.8	0.3	0.005
7.260		0.8	0.3	0.006
7.430		0.9	0.3	0.006
7.508		3.3	1.0	0.022
7.706		14853.2	2838.1	99.886