



Agilent Technologies

Sample name: M4847-Tafamidis-01

Instrument: A010

Inj. volume: 3

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B,Sig=210.0,4.0 Ref=off

RT [min]	Area	Height	Area%
7.617	33.24	9.70	0.27
9.062	42.87	10.14	0.35
9.229	36.33	9.62	0.30
9.442	9.08	1.95	0.08
12.598	11976.41	1371.44	99.00