

Sample Name : BAY 87-2243-01

Method Info : Column:SunFire C18 (4.6x 50 mm,3.5 um)


Mobile phase: H2O(0.01%TFA) (A) / ACN(0.01%TFA) (B)

Elution program: Gradient from 10 to95% of B in 1.4 min at 2ml/min


Temperature: 50 °C

Detection: UV (214 , 4 nm) and MS (ESI, Pos mode ,110 to 1000 amu)


*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.67	0.03	1096.80	472.05	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.67	0.03	2772.04	1154.91	97.72
2.10	0.01	64.75	69.45	2.28

Ret. Time: 1.67

<<<< POSITIVE SPECTRA >>>>



Ret. Time: 2.10

