



Sample name: M4890-JZL 195-01

Instrument: A012

Inj. volume: 1.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B,Sig=210,4 Ref=off

RT [min]	Area	Height	Area%
6.369	2.51	0.93	0.03
7.943	9292.64	1947.26	99.73
9.510	2.17	0.71	0.02
9.615	1.45	0.47	0.02
9.840	18.84	5.02	0.20