



Agilent Technologies

Sample name: M4911-Isoxazole 9-02-1

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1G,Sig=280.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
6.300	BB	0.19	4.04	1.28	0.05
6.886	MM m	0.15	0.57	0.17	0.01
7.352	BB	0.23	7.36	2.12	0.08
7.581	MM m	0.14	0.23	0.05	0.00
7.713	MM m	0.13	1.72	0.55	0.02
8.100	BB	0.67	8795.93	2350.08	99.84