

Sample Name: M4940-SC75741-02

Acq. Operator : Z

Acq. Instrument : instr1


Location : Vial 12

Inj Volume : 3.000 µl

Acq. Method : D:\CHEM32\1\METHODS\Z.M

Analysis Method : D:\CHEM32\1\METHODS\Z2.M

Additional Info : Peak(s) manually integrated



Area Percent Report

Sorted By : Signal
Multiplier: : 1.0000
Dilution: : 1.0000
Use Multiplier & Dilution Factor with ISTDs

Signal 1: VWD1 A, Wavelength=220 nm

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	7.363	BV	0.0631	10.31991	2.45048	0.1432
2	7.529	VB	0.0819	11.54481	2.04486	0.1602
3	9.486	BB	0.1275	7185.38770	765.49335	99.6966

Totals : 7207.25242 769.98868

*** End of Report ***