



Sample Name : M5069-S0859-01
Acq. Operator : C22
Spec. Reported : MS Integration
Acq. Method : D:\METHODS\1-POS-3MIN.M
Analysis Method : D:\METHODS\1-POS-3MIN.M
Sample Info : Easy-Access Method: '1-POS-3MIN.M'
Method Info : Column:Xbridge C18(2) (4.6x 50 mm, 3.5um)
Mobile phase: H2O(10 mmol NH4HCO3) (A) / ACN(B)
Elution program: Gradient from 5 to 95% of B in1.5min at 1.8ml/min
Temperature: 50 °C
Detection: UV (214 , 4 nm) and MS (ESI, Pos mode ,103 to 800 amu)


*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 80



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.43	0.02	6.59	5.56	0.20
1.50	0.02	15.91	14.35	0.47
1.54	0.02	3321.66	2410.79	99.09
2.04	0.02	7.90	5.76	0.24


Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.54	0.02	670.74	509.39	100.00

Ret. Time: 1.55


<<<< POSITIVE SPECTRA >>>>

*ES-API Positive



Ret. Time: 2.04

*ES-API Positive



Ret. Time: 2.29

*ES-API Positive

