



Agilent Technologies

Sample name: M5084-Raltegravir potassium-01

Injection date:

Instrument: A012

Location:

Inj. volume: 2


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1A,Sig=220.0,4.0 Ref=off

| RT [min] | Peak<br>Resolution<br>USP | Area   | Height | Area%  |
|----------|---------------------------|--------|--------|--------|
| 6.980    |                           | 2.1    | 0.6    | 0.026  |
| 7.320    |                           | 4.9    | 1.0    | 0.061  |
| 7.567    |                           | 11.0   | 3.7    | 0.136  |
| 7.697    |                           | 8032.4 | 2153.5 | 99.777 |