



Agilent Technologies

Sample name: M5139-Amlodipine-01

Instrument: A010

Inj. volume: 3.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
6.235	BB	0.12	3.44	1.16	0.04
6.584	MB m	0.27	42.34	8.78	0.46
6.812	BB	0.13	3.83	1.23	0.04
7.377	BM m	0.42	9144.04	2045.36	99.42
7.853	MM m	0.12	0.96	0.23	0.01
9.224	MM m	0.10	1.38	0.41	0.01
9.423	MM m	0.12	1.62	0.44	0.02