



Sample name: M5143-Stenabolic (SR9009)-02

Injection date:

Instrument: 2.5

Location:

Inj. volume: HPLC


Type:

Acq. method:

Sample amount:

Processing method:

Manually modified:



Signal: DAD1A,Sig=220,4 Ref=off

RT [min]	Peak Resolution USP	Width [min]	Area	Height	Area%
9.485		0.06	9.2	2.4	0.103
10.371		0.49	8789.1	1777.7	98.120
10.908		0.07	33.8	8.3	0.377
11.808		0.11	71.5	11.4	0.798
14.231		0.07	53.9	12.2	0.601
Sum			8957.47		