


SAMPLE INFORMATION

Sample Name:	M5174-Icariin-02	Acquired By:	System
Sample Type:	Unknown	Sample Set Name:	
Vial:	14	Acq. Method Set:	Icariin_KHFF
Injection #:	2	Processing Method:	Icariin
Injection Volume:	10.00 ul	Channel Name:	270.0nm
Run Time:	60.0 Minutes	Proc. Chnl. Descr.:	PDA 270.0 nm

Auto-Scaled Chromatogram



Peak Results

	Name	RT	Area	Height	% Area
1		8.950	47714	3087	0.67
2		28.950	7040194	201942	98.29
3		38.543	74772	11238	1.04

Mobile phase A:Acetonitrile

Mobile phase B:Water

Gradient programme(T/%A):0/24,30/26,40/90,45/90,45.1/24,52/24

Flow rate:1.0ml/min

Column temp:30°

Chromatographic column:Kromasil C18 250*4.6 5um