



Sample name: M5246-EW-7197-02

Instrument: A012

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1E,Sig=240,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
4.086	MM m	0.38	28.80	2.77	0.28
4.356	MM m	0.12	0.97	0.34	0.01
4.616	MM m	0.09	1.86	0.67	0.02
5.087	MM m	0.23	3.36	0.70	0.03
5.303	MM m	0.12	1.54	0.50	0.01
5.462	MM m	0.14	1.08	0.25	0.01
5.605	MM m	0.19	1.98	0.41	0.02
5.839	BB	0.16	9.89	3.04	0.10
5.938	BB	0.08	2.37	1.01	0.02
6.063	BB	0.39	10263.83	3013.18	99.03
6.666	MV m	0.19	41.39	13.18	0.40
6.872	BB	0.12	7.44	2.51	0.07