


Sample Name : M5247-BAY-61-3606-01  
Acq. Operator : 9  
Spec. Reported : Integration

Location : P1-F-01  
Inj : 0  
Inj Volume : 6 ul



Integration Results for DAD1 A, Sig=214,4 Ref=off

| RetTim | Width | Area    | Height  | Area% |
|--------|-------|---------|---------|-------|
| 2.37   | 0.04  | 9.26    | 3.35    | 0.14  |
| 2.57   | 0.05  | 7.20    | 2.41    | 0.11  |
| 3.26   | 0.07  | 6376.18 | 1428.59 | 99.25 |
| 3.80   | 0.05  | 31.93   | 9.35    | 0.50  |