



Agilent Technologies

Sample name: M5249-Sotagliflozin-02

Instrument: A012

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
7.312	BB	0.14	4.65	1.12	0.05
7.475	BB	0.16	7.21	1.84	0.08
8.720	MM m	0.49	9435.78	1835.49	99.85
9.659	MM m	0.13	2.54	0.69	0.03