



Agilent Technologies

Sample name: M5326-HTH-01-015-02

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
6.009	BB	0.17	4.93	1.33	0.05
6.132	BB	0.09	17.75	6.68	0.16
6.249	BB	0.32	10801.31	2618.22	98.93
6.579	BB	0.22	58.55	9.69	0.54
6.831	MM m	0.10	5.20	1.42	0.05
6.902	MB m	0.15	14.93	3.60	0.14
7.199	MB m	0.15	6.39	1.66	0.06
8.803	MB m	0.13	3.81	1.05	0.03
9.774	MB m	0.15	5.30	1.26	0.05