

LC/MS Report

Sample Name : M5335-PF-4989216-01

Method Info : Mobile Phase : A:water(0.01%TFA) B:ACN(0.01%TFA)

Gradient:5%B increase to 95%B within 1.5min,95%B for 1.5min


Flow Rate :2.0ml/min

Column :SunFire C18 3.5um, 4.6mm*50mm


Column Temperature:50C

Detection:UV(214nm,4nm) and MS(ESI,POS Mode,110-1300amu)


*DAD1 A, Sig=2144 Ref=off



*DAD1 B, Sig=2544 Ref=off



MSD1 TIC, MS File ES-AP I, Pos, Scan, Frag: 100



Integration Results for DAD1 A, Sig=2144 Ref=off

RetTim	Width	Area	Height	Area%
1.49	0.01	1381.35	1634.24	100.00

Integration Results for DAD1 B, Sig=2544 Ref=off

RetTim	Width	Area	Height	Area%
1.49	0.01	1158.55	1370.54	100.00

LC/MS Report

Ret. Time: 1.50

