

LC-MS Report


Sample Name : M5338-PS-1145-0

Acq Operator : LCMS-13

Inj. Vol. : 5ul

Location : Vial 2

Method Info. : column:C18;column size:4.6*50mm ;mobile phase:B (ACN)
A(0.02%NH4Ac); gradient(B%)!as Acq.Method above



Integration Results

Signal 1: DAD1 A, Sig=214,8 Ref=off


#	R.T.	Type	Height	Height%	Width	Area	Area %
1	3.039	MF	2877.416	100.000	0.056	9743.443	100.000

Signal 2: DAD1 B, Sig=254,16 Ref=off

#	R.T.	Type	Height	Height%	Width	Area	Area %
1	3.039	MF	2940.905	100.000	0.055	9735.793	100.000

Signal 3: MSD1 TIC, MS File

#	R.T.	Type	Height	Height%	Width	Area	Area %
1	3.063	MM	2.373e+006	100.000	0.090	1.283e+007	100.000



*** End of Report ***