


Data file: D:\CHEM32\1\DATA\240810\M54237Final.D
Sample name: M54237-Kisspeptin-10, human-02
Description: Mobile phase: Buffer A: 0.1%TFA in H2O
Buffer B: 0.1%TFA in ACN
Gradient: 27-37% B in 20min, Flow rate: 1ml/min
Column: Unitary C18, 4.6*150mm, 5um 100A
Instrument: Agilent HPLC 1100

Instrument: 1100
Injection date: 2024/8/10 14:31:32
Last changed: 2024/8/10 14:32:23

Location: 51
Injection: 1 of 1
Injection volume: 10.000
Acq.operator: SYSTEM



Signal: DAD1 A, Sig=220,4 Ref=360,100

RT[min]	Type	Width[min]	Area	Height	Area%	Name
7.454	BV E	0.2317	6.7108	0.4457	0.0371	
7.832	VV E	0.2453	66.3677	4.0908	0.3669	
8.435	VV R	0.2687	17940.9473	1063.9946	99.1955	
9.590	VV E	0.2337	24.5799	1.6133	0.1359	
9.908	VV E	0.1313	3.1546	0.3706	0.0174	
10.299	VB E	0.2399	5.7592	0.4201	0.0318	
10.843	BB	0.2238	38.9340	2.7069	0.2153	
Sum			18086.4534			