

MASS SPECTROMETRY REPORT



Sample Description

Analyzed date: 2024-08-10
Analyst: SHAO
Sample: M54237-Kisspeptin-10, human-02
M.W.: 1301.6
Lot. No.:

Instrument

Probe:
Nebulizer Gas Flow:
CDL:
CDL Temp.:
Block Temp.:

Waters 2695

ESI
1.5L/min
-20.0v
250 °C
200 °C

Probe Bias: +4.5kv
Detector: 1.5kv
T. Flow: 0.5ml/min
B. Conc.: 50%H2O/50%ACN