



Agilent Technologies

Sample name: M55326-JG-2016-02

Instrument: A010

Inj. volume: 2.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1A, Sig=220.0, 4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
6.261		4.7	1.4	0.046
6.645		0.7	0.2	0.007
7.108		1.3	0.4	0.013
7.272		3.8	1.2	0.037
7.420		20.6	5.9	0.203
7.635		1.1	0.4	0.011
7.846		2.6	0.7	0.026
7.926		4.6	1.0	0.045
9.145		10059.6	1929.4	98.758
9.264		79.4	22.9	0.780
9.967		6.6	1.9	0.065
10.479		1.1	0.3	0.011