



8.310
8.305
8.292
8.287
8.282
8.007
8.001
7.996
7.983
7.979
7.973
7.467
7.421
7.412
6.843
6.834

3.725

3.342

2.518

2.514

2.510

2.505

2.501

Current Data Parameters
NAME Dantrolene sodium hemiheptahydrate-01
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters

Time 10.19
INSTRUM spect
PROBHD 5 mm PADUL 13C
PULPROG zg30
TD 65536
SOLVENT DMSO
NS 16
DS 2
SWH 8012.820 Hz
FIDRES 0.122266 Hz
AQ 4.0894465 sec
RG 199.46
DW 62.400 usec
DE 6.50 usec
TE 297.0 K
D1 1.00000000 sec
TD0 1

===== CHANNEL f1 =====
SFO1 400.1324710 MHz
NUC1 1H
P1 10.00 usec
PLW1 20.00000000 W

F2 - Processing parameters
SI 65536
SF 400.1300000 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00

9.5 9.0 8.5 8.0 7.5 7.0 6.5 6.0 5.5 5.0 4.5 4.0 3.5 3.0 2.5 2.0 1.5 1.0 ppm

2.07
2.10
1.04
1.03
1.00

2.15