



Sample name: M5588-Docetaxel Trihydrate-01

Injection date:

Instrument: A012

Location:

Inj. volume: 5


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1D,Sig=230,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
6.514	MM m	0.09	2.15	0.75	0.02
8.173	MV m	0.14	5.25	1.34	0.06
8.962	BB	0.48	8973.22	1703.77	99.15
9.524	MM m	0.12	1.92	0.55	0.02
9.882	MB m	0.34	63.64	11.78	0.70
10.287	MM m	0.15	3.76	0.92	0.04