


HPLC REPORT

Sample: M58139-CP05-02 Analyzed date: 2024-10-31
Analyst: Apeptide
Column: 4.6×250mm,Sinobrom ODS-BP 5μm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient: A B
0.0min 10% 90%
25.0min 35% 65%
25.1min 100% 0%
30.0min Stop
Volume: 5μl
Wavelength: 220nm
Flow rate: 1.0ml/min



Peak	Time	Height	Area	Conc.
1	12.510	316.113	1753.496	0.0483
2	13.100	13020.482	111956.641	3.0810
3	13.403	217176.141	3481484.750	95.8080
4	17.145	1055.552	8032.086	0.2210
5	17.810	1674.125	13126.112	0.3612
6	24.007	2223.881	17461.004	0.4805

Total 100.000