



Agilent Technologies

Sample name: M5878-Pioglitazone HCl-02

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
5.951	MM m	0.15	1.46	0.42	0.02
6.275	BM m	0.29	8814.26	2481.53	99.86
6.503	MM m	0.09	1.77	0.41	0.02
6.694	BB	0.10	5.13	1.93	0.06
6.786	BB	0.10	4.14	1.46	0.05