



Agilent Technologies

Sample name: M59015-(R)-Zunsemetinib-02

Injection date:

Instrument: A012

Location:

Inj. volume: 4


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1A, Sig=220,4 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
5.639		1.3	0.5	0.010
5.811		0.6	0.3	0.005
5.889		3.8	1.0	0.030
6.051		1.2	0.4	0.009
6.146		4.7	1.7	0.036
6.345		19.3	6.2	0.151
6.878		1.9	0.7	0.015
7.077		4.1	1.7	0.032
7.152		7.6	3.0	0.059
7.275		1.7	0.6	0.013
7.435		12674.3	2593.6	98.836
7.645		20.6	6.1	0.161
7.923		78.6	20.5	0.613
9.578		0.0	0.0	0.000
9.873		3.8	1.0	0.030