




Sample name: M6004-Testosterone propionate-03 Injection date:
Instrument: 4 Location:
Inj. volume: HPLC Type:
Acq. method: Sample amount:
Processing method:
Manually modified:



Signal: DAD1E,Sig=240,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%	Name
8.037	MM m	0.06	0.57	0.16	0.00	
8.371	MM m	0.05	0.33	0.12	0.00	
8.856	MM m	0.05	0.88	0.28	0.00	
11.037	MM m	0.06	2.11	0.56	0.01	
13.518	BB	1.33	27233.71	2549.19	99.99	
	Sum		27237.60			