

# Shimadzu LCMS-029 Report

## Sample Information

### LCMS


Acquired by : System Administrator  
Sample Name : M6143-SP2509-01  
Method File : B 5%-95%-P-new - 2.6min.lcm

Date Acquired : 2021/8/21 10:14:26  
Injection Volume : 2  
Vial : 74


Column : Waters SunFire C18 50\*4.6mm 5um 2.000ml/min 2.6min  
Pump A : 0.1% FA and 10% ACN in H2O  
Pump B : 0.1% FA and 10% H2O in ACN

## Chromatogram


mAU



mAU



段号1



## Integration Results for DAD

### PDA Ch1 254nm

| Peak# | Ret. Time | Height | Area   | Area%   |
|-------|-----------|--------|--------|---------|
| 1     | 1.801     | 556770 | 966437 | 100.000 |
| 总计    |           | 556770 | 966437 | 100.000 |


### PDA Ch2 214nm

| Peak# | Ret. Time | Height  | Area    | Area%   |
|-------|-----------|---------|---------|---------|
| 1     | 1.801     | 1703070 | 2901972 | 100.000 |


| Peak# | Ret. Time | Height  | Area    | Area%   |
|-------|-----------|---------|---------|---------|
| 总计    |           | 1703070 | 2901972 | 100.000 |

## Mass Spectrum


Rt=0.754min 正




Rt=1.384min 正




Rt=1.547min 负




Rt=1.542min 正




Rt=1.750min 负




Rt=1.767min 正




Rt=2.110min 负



Rt=2.217min 正



Rt=2.245min 负



Rt=2.402min 负

