




Sample Name : M6151-LOXO-101-07
Acq. Operator : C2
Spec. Reported : MS Integration
Acq. Method : D:\METHODS\2-POS-MON-1.M
Analysis Method : D:\METHODS\2-POS-MON-1.M
Sample Info : Easy-Access Method: '2-POS-MON-1'
Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)
 Gradient: 5%-95% B in 1.5min,95%B for 1.5 min
 Flow Rate:2ml/min
 Column:Sunfire ,4.6*50mm,3.5um
 Oven Temperature: 50 C
 MASS Range:120-1000



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.62	0.01	981.78	1131.65	99.12
2.55	0.02	8.76	7.58	0.88

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.62	0.01	1170.85	1349.83	100.00

Ret. Time: 1.64

<<<< POSITIVE SPECTRA >>>>

ES-API Positive

