




Sample name: M6178-Tucidinostat (Chidamide)-02 **Injection date:**
Instrument: A010 **Location:**
Inj. volume: 2 **Type:**
Acq. method: HPLC **Sample amount:**
Processing method:
Manually modified:



Signal: DAD1C,Sig=254.0,4.0 Ref=off

RT [min]	Peak Resolution USP	Width [min]	Area	Height	Area%
4.959		0.16	171.5	55.7	1.85
5.474		0.40	9074.3	2699.2	98.06
5.820		0.12	7.6	2.0	0.08