



Agilent Technologies

Sample name: M6228-Wilforine-02

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
7.699	MM m	0.08	3.92	1.30	0.03
7.770	MM m	0.08	12.50	3.36	0.09
7.854	MM m	0.12	27.26	8.04	0.20
7.991	BB	0.13	7.75	2.53	0.06
8.394	BM m	0.58	13772.27	2595.68	99.46
8.638	MB m	0.18	11.07	2.52	0.08
10.876	MB m	0.16	12.37	2.88	0.09