




LCMS- Report

LCMS
Acquired by : SD
Sample Name : M6260-ROC-325-02
Method File : 3.lcm

Method info. : C18 4.6*50mm 5um 2.000ml/min 2.6min Column Temperature: 40 °C
Gradient: 5% B hold for 0.2min, increase to 95 % B within 1.40 min, hold at 95 % B for 0.9 min, then back to 5% B within 0.01 min
Pump A : 0.03% TFA in H2O
Pump B : 0.03% TFA in ACN

Chromatogram



Integration Results for DAD

PDA Ch1 254nm

Peak#	Ret. Time	Height	Area	Area%
1	1.383	809137	3775204	100.000
Total		809137	3775204	100.000

PDA Ch2 214nm

Peak#	Ret. Time	Height	Area	Area%
1	1.383	716791	2965609	100.000
Total		716791	2965609	100.000

Mass Spectrum

Rt=1.387min Positive

