



Sample name: M6297-Upadacitinib-01

Injection date:

Instrument: A010

Location:

Inj. volume: 3


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
4.456	MM m	0.09	0.57	0.18	0.01
4.540	MM m	0.10	1.86	0.65	0.03
4.663	MM m	0.08	2.65	0.87	0.04
4.737	MM m	0.09	2.08	0.62	0.03
4.887	MM m	0.10	1.45	0.51	0.02
5.277	MM m	0.12	4.46	1.49	0.06
5.409	MM m	0.10	0.49	0.18	0.01
5.579	BB	0.17	10.82	2.41	0.16
5.740	BB	0.35	6926.78	1645.40	99.60
6.289	MM m	0.09	1.75	0.66	0.03
6.441	MM m	0.07	0.97	0.39	0.01
6.920	MM m	0.08	0.85	0.32	0.01