



Agilent Technologies

Sample name: M6499-Bepiridil hydrochloride-02

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B, Sig=210.0, 4.0 Ref=off

RT [min]	Area	Height	Area%
7.365	27.17	8.68	0.29
7.762	3.24	0.60	0.03
9.292	9299.62	1267.83	99.64
9.458	3.61	1.84	0.04