



Agilent Technologies

Sample name: M6499-Bepidil hydrochloride-03

Injection date:

Instrument: A012

Location:

Inj. volume: 2


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1H,Sig=200.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
9.407	BM m	0.49	9289.92	1160.89	99.87
10.160	MM m	0.11	11.67	3.30	0.13