




Sample Name : M6852-JTE 607 dihydrochloride-02
Acq. Operator : 11
Spec. Reported : MS Integration
Acq. Method : D:\METHODS\1-POS-3MIN.M
Analysis Method : D:\METHODS\1-POS-3MIN.M
Sample Info : Easy-Access Method: '1-POS-3MIN.M'
Method Info : Column:Xbridge C18(2) (4.6x 50 mm, 3.5um)
 Mobile phase: H2O(10 mmol NH4HCO3) (A) / ACN(B)
 Elution program: Gradient from 5 to 95% of B in1.5min at 1.8ml/min
 Temperature: 50 °C
 Detection: UV (214 , 4 nm) and MS (ESI, Pos mode ,103 to 800 amu)



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.29	0.02	2056.84	1590.08	98.53
1.41	0.02	30.77	23.77	1.47


Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.29	0.02	685.87	530.82	98.78
1.41	0.02	8.47	6.56	1.22

Ret. Time: 1.31

<<< POSITIVE SPECTRA >>>

*ES-API Positive



Ret. Time: 2.03

*ES-API Positive

