



Agilent Technologies

Sample name: M7189-Resiquimod-02

Instrument: A010

Inj. volume: 1

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
4.951	MM m	0.10	1.26	0.47	0.01
5.134	BB	0.18	15.69	5.22	0.18
5.326	BB	0.15	16.34	5.29	0.19
5.439	MM m	0.11	0.66	0.27	0.01
5.538	MM m	0.09	0.76	0.26	0.01
5.854	BB	0.39	8710.30	2559.20	99.28
6.208	BB	0.14	10.22	2.75	0.12
6.519	BB	0.10	6.27	2.25	0.07
6.741	BB	0.14	3.99	1.05	0.05
6.976	MM m	0.11	1.98	0.59	0.02
7.476	MM m	0.20	3.77	0.60	0.04
9.408	MM m	0.12	2.40	0.65	0.03