




Sample Name : M7353-Tamoxifen-09
Acq. Operator : 11
Spec. Reported : MS Integration
Acq. Method : D:\METHODS\2-POS-MON-1.M
Analysis Method : D:\METHODS\2-POS-MON-1.M
Sample Info : Easy-Access Method: '2-POS-MON-1'
Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)
Gradient: 5%-95% B in 1.5min, 95%B for 1.5 min
Flow Rate: 2ml/min
Column: Sunfire, 4.6*50mm, 3.5um
Oven Temperature: 50 C
MASS Range: 120-1000



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.70	0.01	6.95	9.08	0.17
1.75	0.02	4146.43	2799.94	98.91
2.52	0.02	9.21	8.54	0.22
2.84	0.03	29.33	12.39	0.70

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.75	0.02	3076.86	2797.79	100.00

Ret. Time: 1.77 <<<< POSITIVE SPECTRA >>>>

