



Agilent Technologies

Sample name: M7353-Tamoxifen-12

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B,Sig=210.0,4.0 Ref=off

RT [min]	Area	Height	Area%
9.450	7.09	2.40	0.03
9.879	20596.08	2017.64	99.78
13.852	16.40	3.33	0.08
16.519	21.19	4.81	0.10